Design Guidelines for Hospitals and Day Procedure Centres

Standard Components Room Data Sheets

Room Name

Catheter Lab Control/ Reporting Room

Room Code CORT

Area M2

10

Description Special Requirement

Room adjacent to Catheter Laboratory room/s with direct visibility to the patient for remote control of equipment, and review and reporting of procedure images. The Control/ Reporting room may serve two procedure rooms.

Acoustics: Acoustic privacy is required

Additional Design Considerations:

- Radiation shielding as advised by Radiation Consultant
- Cable ducts from Catheter Laboratory and Computer/Equipment room may be required

Room	Fabric and out	line specific	ations		Remarks
1000	Floor finish	Vinyl	Standard slip resistant	Seamless, coved	
2000	Skirting	Vinyl	Prefinished	Floor vinyl coved, 150mmH	
3000	Wall finish	Paint	Acrylic, washable		
5060	Ceiling	Acoustic	Prefinished	Drop-in tiles, 600x1200	
6030	Cornice	Aluminium	Powdercoat	24mm shadow angle	
7025	Observation	Lead glass	Clear	Internal window	to procedure room
8000	Door protection	Vinyl	Prefinished	Kickplate, 300mm high	
9004	Doors	Solid Core/ Glass	Paint	Single leaf, observation panel	staff entry/ exit

Fittings and Furniture

Fitting ID	Category	Fitting	Grp	Quantity	Selection
1050	Furniture/ Fitting	Bench: laminate	1	1	
3400	Metalwork	Coat hook	1	2	
17360	Furniture/ Fitting	Pinboard, fabric covered	1	1	over bench
20000	Furniture/ Fitting	Shelving unit, adjustable	1	1	for books, manuals, over part bench

Fixtures, Equipment and associated services

FFE ID	Description	Grp	Quantit	Ele	CdW	HtW	WmW	Тар	Dns	Gas	Stm	Selection
62150	XR control console	2	1	yes								Angiography equipment

Services

Service ID	Service Category	Description	Quantity	Remarks
1006	Communications	Voice / Data outlet - double	4	
1040	Communications	Intercom	1	with handsfree microphone to procedure room
5000	HVAC	Heating, Ventilation & Airconditioning		Refer to Part E Building Services
6000	Lighting	General		
6055	Lighting	Special: task light- built in		under O/H shelving
9001	Power	GPO: Double	9	
9060	Power	Special: for equipment	1	Angiography equipment





